

# ALL CUP TESTS

## WITH CUT-OFF LEVELS

|                                 | ID-CP6    | ID-US6 | ID-US6-BUP | ID-H7-1 | ID-H9-1 | ID-CP10 | ID-CP12-BUP | ID-CP12-BUP-ADULT | ID-US12-BUP | ID-H12-EFKK-1 | ID-H12-1 | ID-H12-2 | ID-H13-1 | ID-US14-ADULT | ID-H14-ADULT-EFKT-2 | ID-H16-ADULT-EFKT | ID-US16-ADULT-EFKT |
|---------------------------------|-----------|--------|------------|---------|---------|---------|-------------|-------------------|-------------|---------------|----------|----------|----------|---------------|---------------------|-------------------|--------------------|
| SUBSTANCE TESTED                | CUP TESTS |        |            |         |         |         |             |                   |             |               |          |          |          |               |                     |                   |                    |
| Amphetamines (AMP)              | 1000      | 500    |            | 500     | 500     | 1000    | 1000        | 1000              | 500         | 500           | 500      | 500      | 500      | 500           | 500                 | 500               | 500                |
| Barbiturates (BAR)              |           |        |            |         | 300     | 300     | 300         | 300               | 300         |               | 300      | 300      | 300      | 300           |                     | 300               | 300                |
| Benzodiazepines (BZO)           | 300       | 300    | 300        | 300     | 300     | 300     | 300         | 300               | 300         | 300           | 300      | 300      | 300      | 300           | 300                 | 300               | 300                |
| Buprenorphine (BUP)             |           |        | 10         | 10      | 10      |         | 10          | 10                | 10          | 10            | 10       | 10       | 10       | 10            | 10                  | 10                | 10                 |
| Cocaine (COC)                   | 300       | 150    | 150        | 150     | 150     | 300     | 300         | 300               | 150         | 150           | 150      | 150      | 150      | 150           | 150                 | 150               | 150                |
| Ecstasy (MDMA)                  |           |        |            |         |         |         | 500         | 500               | 500         |               | 500      | 500      | 500      | 500           | 500                 | 500               | 500                |
| Marijuana (THC)                 | 50        | 50     |            | 50      |         | 50      | 50          | 50                | 50          | 50            | 50       | 50       | 50       | 50            | 50                  | 50                | 50                 |
| Methadone (MTD)                 |           |        | 300        |         | 300     | 300     | 300         | 300               | 300         |               | 300      | 300      | 300      | 300           | 300                 |                   |                    |
| Methadone Metabolite (EDDP)     |           |        |            |         |         |         |             |                   |             |               |          |          |          | 300           |                     | 300               | 300                |
| Methamphetamines (MET)          |           |        |            |         | 500     | 1000    | 1000        | 1000              | 500         | 500           | 500      | 500      | 500      | 500           | 500                 | 500               | 500                |
| Opiates (MOP) or (OPI300)       | 300       | 300    | 300        | 300     | 300     | 300     | 300         | 300               | 300         | 300           | 300      | 300      | 300      | 300           | 300                 | 300               | 300                |
| Oxycodone (OXY)                 | 100       | 100    | 100        | 100     | 100     | 100     | 100         | 100               | 100         | 100           | 100      | 100      | 100      | 100           | 100                 | 100               | 100                |
| Phencyclidine (PCP)             |           |        |            |         |         | 25      | 25          | 25                | 25          |               | 25       |          | 25       | 25            |                     | 25                | 25                 |
| Tricyclic Antidepressants (TCA) |           |        |            |         |         |         |             |                   |             |               |          | 1000     | 1000     | 1000          |                     |                   |                    |
| ** Cotinine / Nicotine (COT)    |           |        |            |         |         |         |             |                   |             |               |          |          |          |               |                     |                   |                    |
| ** Ethyl Glucuronide (ETG)      |           |        |            |         |         |         |             |                   | 500         |               |          |          |          |               | 500                 | 500               | 500                |
| ** Fentanyl (FEN)               |           |        |            |         |         |         |             |                   | 20          |               |          |          |          |               | 20                  | 25                | 25                 |
| ** Gabapentin (GAB)             |           |        |            |         |         |         |             |                   |             |               |          |          |          |               |                     |                   |                    |
| ** Kratom (KRA)                 |           |        |            |         |         |         |             |                   | 100         |               |          |          |          |               |                     |                   |                    |
| ** Synthetic Marijuana (K2)     |           |        |            |         |         |         |             |                   | 25          |               |          |          |          |               | 25                  | 25                | 25                 |
| ** Tramadol (TRA)               |           |        |            |         |         |         |             |                   |             |               |          |          |          |               | 100                 | 100               | 100                |
| *ADULTERANTS                    |           |        |            |         |         |         |             |                   |             |               |          |          |          |               |                     |                   |                    |
| Creatinine (C)                  |           |        |            |         |         |         |             | ✓                 |             |               |          |          |          | ✓             | ✓                   | ✓                 | ✓                  |
| Glutaraldehyde (G)              |           |        |            |         |         |         |             |                   |             |               |          |          |          | ✓             |                     | ✓                 | ✓                  |
| Acidic or Alkaline pH (P)       |           |        |            |         |         |         |             | ✓                 |             |               |          |          |          | ✓             | ✓                   | ✓                 | ✓                  |
| Nitrite (N)                     |           |        |            |         |         |         |             |                   |             |               |          |          |          | ✓             |                     | ✓                 | ✓                  |
| Oxidants (O)                    |           |        |            |         |         |         |             |                   |             |               |          |          |          | ✓             | ✓                   | ✓                 | ✓                  |
| Specific Gravity (S)            |           |        |            |         |         |         |             | ✓                 |             |               |          |          |          | ✓             |                     | ✓                 | ✓                  |

\*\* NOTE: Premium panels not CLIA Waived intended for forensic use only. \* NOTE: Urine tests with adulterants measure levels of substances in human urine which may indicate sample tampering.

A green box indicates the substance or adulterant in the far left column is tested for.

The number inside the box indicates the cut-off level in ng/ml. A cut-off level is the concentration a med or narcotic must exceed in the urine sample to trigger a positive result.

Tests above with the American Flag are designed & manufactured in the United States.

All Identify Diagnostics & Identify Diagnostics USA tests are made with USA reagents.

Reagents are the substances on each panel strip which produce the urine test result.

# ALL DIP TESTS

## WITH CUT-OFF LEVELS

| SUBSTANCE TESTED                | DIP TESTS      |                |                |                |                |                |            |            |             |               |
|---------------------------------|----------------|----------------|----------------|----------------|----------------|----------------|------------|------------|-------------|---------------|
|                                 | ID-US1-DIP-COT | ID-US1-DIP-THC | ID-US1-DIP-ETG | ID-US1-DIP-FEN | ID-US1-DIP-GAB | ID-US1-DIP-KRA | ID-CP5-DIP | ID-CP6-DIP | ID-CP10-DIP | ID-US10-DIP-1 |
| Amphetamines (AMP)              |                |                |                |                |                | 1000           |            | 1000       | 500         | 1000          |
| Barbiturates (BAR)              |                |                |                |                |                |                |            | 300        | 300         | 300           |
| Benzodiazepines (BZO)           |                |                |                |                |                |                | 300        | 300        | 300         | 300           |
| Buprenorphine (BUP)             |                |                |                |                |                |                |            |            | 10          | 10            |
| Cocaine (COC)                   |                |                |                |                |                | 300            | 300        | 300        | 150         | 300           |
| Ecstasy (MDMA)                  |                |                |                |                |                |                |            | 500        |             | 500           |
| Marijuana (THC)                 | 50             |                |                |                |                | 50             | 50         | 50         | 50          | 50            |
| Methadone (MTD)                 |                |                |                |                |                |                |            | 300        | 300         | 300           |
| Methadone Metabolite (EDDP)     |                |                |                |                |                |                |            |            |             |               |
| Methamphetamines (MET)          |                |                |                |                |                | 1000           | 1000       | 1000       | 500         | 1000          |
| Opiates (MOP) or (OPI300)       |                |                |                |                |                | 300            | 300        | 300        | 300         | 300           |
| Oxycodone (OXY)                 |                |                |                |                |                |                | 100        | 100        | 100         | 100           |
| Phencyclidine (PCP)             |                |                |                |                |                |                |            |            |             | 25            |
| Tricyclic Antidepressants (TCA) |                |                |                |                |                |                |            |            |             |               |
| ** Cotinine / Nicotine (COT)    | 200            |                |                |                |                |                |            |            |             |               |
| ** Ethyl Glucuronide (ETG)      |                | 500            |                |                |                |                |            |            |             |               |
| ** Fentanyl (FEN)               |                |                | 25             |                |                |                |            |            |             |               |
| ** Gabapentin (GAB)             |                |                |                | 2000           |                |                |            |            |             |               |
| ** Kratom (KRA)                 |                |                |                |                | 100            |                |            |            |             |               |
| ** Synthetic Marijuana (K2)     |                |                |                |                |                |                |            |            |             |               |
| ** Tramadol (TRA)               |                |                |                |                |                |                |            |            |             |               |

\*\* NOTE: Premium panel not CLIA Waived intended for forensic use only.

A green box  indicates the substance in the far left column is tested.

The number inside the box indicates the cut-off level in ng/ml. A cut-off level is the concentration a med or substance must exceed in the urine sample to trigger a positive result.

Learn more about cut-off levels & detection times in this image gallery.

Tests above with the American Flag  are designed & manufactured in the United States.

All Identify Diagnostics & Identify Diagnostics USA tests are made with USA reagents.

Reagents are the substances on each panel strip which produce the urine test result.

**Identify®**  
Diagnostics

ID-CP12-BUP



**Identify®**  
Diagnostics USA

ID-US12-BUP




**Identify®**  
Health


ID-H12-1



**5 PANEL**  
(single sided)



**12 PANEL**  
(double sided)



**1 PANEL**  
(FEN)  
fentanyl



**1 PANEL**  
(ETG)  
alcohol



**1 PANEL**  
(COT)  
nicotine



**10 PANEL**  
(double sided)



**Identify®**  
Diagnostics

**Identify®**  
Diagnostics USA

### Target Markets

In Office Physician Labs- CLIA  
Employment Testing  
Parental / Family Testing  
Government  
Addiction

### Distributed By

Medical Distribution Group, Inc.  
2361 Whitfield Park Avenue  
Sarasota, FL 34243  
P: 727-744-2967  
F: 727-499-7355