

Identify[®] CUP TESTS

* Forensic Use Only Panel or Device.

ID-CP5 ID-CP6 ID-H6-1 ID-H7-1 ID-H9-1 ID-CP10 ID-CP12-BUP ID-CP12-BUP-ADULT ID-US12-BUP MADE IN USA ID-H12-EFKK-1 * ID-H12-1 ID-H12-2 ID-H13-1 ID-US14-ADULT MADE IN USA ID-H14-ADULT-1 ID-H14-ADULT-EFKT-2 * ID-H16-ADULT-EFKT *

SUBSTANCE TESTED	CUP TESTS																
Amphetamines (AMP)	1000	500	500	500	1000	1000	1000	500	500	500	500	500	500	500	500	500	500
Barbiturates (BAR)				300	300	300	300	300		300	300	300	300	300		300	
Benzodiazepines (BZO)	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300
Buprenorphine (BUP)			10	10		10	10	10	10	10	10	10	10	10	10	10	10
Cocaine (COC)	300	300	150	150	150	300	300	300	150	150	150	150	150	150	150	150	150
Ecstasy (MDMA)						500	500	500		500	500	500	500	500	500	500	500
Marijuana (THC)	50	50	50	50		50	50	50	50	50	50	50	50	50	50	50	50
Methadone (MTD)				300	300	300	300	300		300	300	300	300	300	300	300	
Methadone Metabolite (EDDP)													300	300		300	
Methamphetamines (MET)				500	1000	1000	1000	500	500	500	500	500	500	500	500	500	500
Opiates (MOP) or (OPI300)	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300	300
Oxycodone (OXY)	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100	100
Phencyclidine (PCP)					25	25	25	25		25		25	25	25		25	
Tricyclic Antidepressants (TCA)											1000	1000	1000	1000			
* Ethyl Glucuronide (ETG)									500							500	500
* Fentanyl (FEN)									25							20	25
* Kratom (KRA)									100								
* Synthetic Marijuana (K2)									25							25	25
* Tramadol (TRA)																100	100

ADULTERANTS

Creatinine (C)									✓						✓	✓	✓	✓
Glutaraldehyde (G)															✓	✓		✓
Acidic or Alkaline pH (P)									✓						✓	✓	✓	✓
Nitrite (N)															✓	✓		✓
Oxidants (O)															✓	✓	✓	✓
Specific Gravity (S)									✓						✓	✓		✓

12 PANEL



Identify[®]
Diagnostics

14 PANEL



Identify[®]
Diagnostics USA

14 PANEL



Identify[®]
Health

Identify[®] DIP TESTS

* Forensic Use Only Panel or Device.

ID-H1-COT-DIP *
 ID-H1-ETG-DIP *
 ID-H1-FEN-DIP *
 ID-H1-GAB-DIP *
 ID-H1-KRA-DIP *
 ID-CP5-DIP
 ID-CP10-DIP
 ID-H10-DIP
 ID-CP12-DIP
 ID-H12-EFT-DIP *

SUBSTANCE TESTED	DIP TESTS								
Amphetamines (AMP)					1000	1000	500	1000	500
Barbiturates (BAR)						300	300	300	
Benzodiazepines (BZO)						300	300	300	300
Buprenorphine (BUP)							10	10	10
Cocaine (COC)						300	300	150	300
Ecstasy (MDMA)							500		500
Marijuana (THC)						50	50	50	50
Methadone (MTD)							300	300	300
Methamphetamines (MET)						1000	1000	500	1000
Opiates (MOP) or (OPI300)						300	300	300	300
Oxycodone (OXY)							100	100	100
Phencyclidine (PCP)									25
* Cotinine / Nicotine (COT)	200								
* Ethyl Glucuronide (ETG)		500							500
* Fentanyl (FEN)			25						50
* Gabapentin (GAB)				2000					
* Kratom (KRA)					100				
* Tramadol (TRA)									100

