



## Product description

Order number	11.2252.001
Product description	V-Monovette® Urine, 4 ml, Vacuum system for closed urine transfer, 75x13mm, round base, paper label, 50/inner box sterile

## Product characteristics

Type of collection	Vacuum
Type of preparation	without preparation
Label/ Print	with paper label
Graduation	Fill mark
Cap	membrane cap
Type of cap	Push cap
Light protection	no
Application	Urine tube
Base shape	round base
Volume of work	4 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	V-Monovette®

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Size**

Sample volume	4 ml
Diameter	13 mm
Length including cap	81 mm
Length excluding cap	75 mm

### **Material & colours**

Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Membrane material	Rubber
Tube material	Polyethylene therephthalate (PET)

### **Purity & certification**

Product in accordance with MDD or IVDD	IVD
CE certified	CE
Purity standard	sterile
Sterilisation	yes
Pyrogen-free	no
Lot or serial number	Lot number

### **Packaging**

Minimum order qty.	500
Pcs. / inner box	50
Pcs. / box	500
Pcs. / pallet	39000
Depth of case	398 mm
Width of case	368 mm
Height of case	130 mm
Case volume	0,019 cbm
Weight of case	3,46 kg
EAN of inner box	4038917244824
EAN of outer box	4038917431842

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

